



## Enjoy a World of Comfort with PlasmaEYEZ

**Now available on all Boston® materials**

- Dramatically reduces the wetting angle
- Removes remaining residues on lens

### Potential Benefits

- Improved comfort
- Enhanced on-eye wettability

### **What is the Plasma Treatment Process?**

Oxygen plasma treatment is a safe and environmentally friendly alternative to traditional cleaning methods. Cold gas plasma is an energetic process utilizing electrical energy to transform the molecular structure of a contact lens surface. This creates desirable surface properties engineered to enhance wettability and comfort.

Products such as contact lenses and intraocular lenses have been successfully modified by plasma treatment to impart protein and cell-repelling characteristics, decrease bacterial adhesion, potentially improve wettability and enhance patient comfort.

**Ask for PlasmaEYEZ on your next GP  
lens order manufactured with Boston materials!**



California  
800-342-2724

Florida  
800-432-3838

Georgia  
800-241-9312

Maryland  
800-221-9235

Minnesota  
800-926-6822

Pennsylvania  
800-245-0797

Washington  
800-426-6241



Minnesota  
800-433-4885



Ohio  
877-633-0300

Iowa  
800-362-1820

Wisconsin  
800-937-4448

Minnesota  
800-736-9339